

## DRK113 Crush Tester



Crush Tester is used to do the Ring Crush test, Edge Crush Test, Flat Press Test, Adhesive Strength Test, and Paper Tube Compress Test.

We also have the accessories for its each test.

Conforms to standards: ISO3035; ISO3037

### **Characters:**

1. Microprocessor control, automatically carry the test according to testing requests.
2. Protect press sensor through the protection switch for over- loading.
3. Analyzing the group of testing data and to show the maximum, minimum, average,
4. standard deviation and coefficient of variation by LCD, print out the data.

### **Technical parameter:**

1. Press measuring range: 60-5000N,
2. Resolution: 1N.
3. Accuracy:  $\pm 1\%$
4. Speed:  $12.5 \pm 2.5$ mm/min
5. Platen distance from top to bottom: 0-150mm (customized for special dimension)
6. Printer: mini-type printer
7. Power : AC220V 50Hz 1A
8. Working condition: Temperature 0-40°C, relative humidity  $< 85\%$
9. Dimension: LxHxW=500x560x330
10. Weight: 70kg
11. Test size