

AA04 – FIRE REACTION WITH HORIZONTAL SAMPLE

TECHNICAL DATASHEET

Model 10091104



Features:

- Stainless steel AISI 304 cabinet
- Cabinet bottom with ventilation slits
- Safety glass door
- 9.5 mm dia. gas burner with support
- Nozzles for methane gas and GPL
- Gas regulation valve
- Extractable horizontal sample support with sample supporting metal wires
- Collecting basin for residual products

General:

- H 360 mm - W 385 mm - D 204 mm.
- Weight 12 kg/13 kg.
- Methane or Propane gas supply

Reference standards

FMVSS	302		
ISO	3795		
DIN	75200		
FAR	Part 25	F part 1	
ABD 0031	AITM 2.003		
BSS	7230		
FIAT	50433		
TL	1010		

Model 10091112



Same features of the standard model, plus:

- **AUTOMATIC timer** for the control of the flame application times



- Facilitates reporting of After-flame values

Code	Description	
10091104	Horizontal test	Standard
10091112	Horizontal test	With built-in timer

